

Paper Chromatography

PURPOSE:

1. To understand and perform a chromatography.
2. To separate the components of commercial food dyes by paper chromatography.

PRINCIPLES:

Chromatography is a method of separating the components of a mixture by distributing them between two phases, one of which is moving past the other one, which is stationary. Since the extent to which the components of the mixture are attracted to the moving phase and respectively to the stationary phase differs, the components of the mixture will be selectively picked up by the moving phase.

The height to which the components of the mixture rise depends upon the relative strength of the:

- mixture and moving phase attractions, and
- mixture and stationary phase attractions.

Thus, the components can be separated in this experiment:

- The **stationary phase** is a good quality **filter paper**.
- The **moving phase** (referred to as the **solvent**) is a 0.1% solution of sodium chloride in water.

The moving phase (solvent) moves along the stationary phase (filter paper) by capillary action due to the attraction between the solvent and the cellulose of which the paper is made.

The mixture (food dye) to be separated is placed on the paper as a small spot, and is carried along with the solvent. Since the different components of the mixture are not attracted to the paper and to the solvent to the same degree, they will move with different speeds along the paper and thus can be separated.


You will analyze three different colors (red, blue and green) of a brand of commercial food dye using a 0.1 % solution of sodium chloride as the solvent. The analysis of the chromatograms will enable you:

- To deduce the extent of relative attractions of the components of the food dyes to the paper
- To deduce the extent of relative attractions of the components of the food dyes to the solvent.

PROCEDURE:


1. Preparation of the filter paper

- The filter paper should be cut into a rectangle about 11 cm x 20 cm.
- Draw a pencil line 1.5 cm from one edge (use a long edge) and another line exactly 8.0 cm above the first. The lower line is your origin or baseline, and the upper line is the desired level for your solvent to reach.
- Place five evenly spaced pencil "X" marks along the lower line, 3 cm from each side and 3.5 cm away from each other.
- Using a pencil, label your chromatogram in the upper right hand side corner with your name and the brand of food dye you are analyzing.
- Using a pencil, label the "X" marks on the baseline with R for Red, B for Blue, G for Green and Numbers for the two unknowns.



2. Preparation of the chromatographic developing chamber.
 - Measure out 30 ml of solvent (0.1% solution of sodium chloride in water) into a 600 ml beaker.
 - Cover the beaker with a watch glass to slow down the evaporation of the solvent.
3. Spotting the filter paper.
 - Use the 3 colors (Red, Blue, and Green) of food dye for analysis.
 - Place 1 drop of each of the followings in three different depressions of your spot plate: Red Food dye, Blue Food dye and Green Food dye.
 - Use an open ended capillary tube to practice applying small spots on a piece of paper towel.
NOTE: Place all used capillary tubes in the common container filled with water provided by your instructor. Try not to break them as they can be used again after proper washing.
 - When you are confident that you master the technique of spotting, place a small spot of each of the three dyes on the X marks.
 - Make sure that the spots are applied exactly over the “X” marks, as to avoid having them below the solvent level when placing the paper in the developing chamber. Use a different capillary tube for each dye and try to keep the spots small (they should be no more than 3 to 5 mm in diameter). If capillary tubing of about 1 mm ID (inside diameter) is used, one application is sufficient. If narrower capillary tubing is used, about 0.4 mm ID, three applications are needed.
 - Let the spots dry between applications to keep the spot size down. Let the paper dry for 5 minutes.

4. Roll the paper into a cylinder and staple the sides without overlap. Place the lower staple above the baseline.
5. Place the cylinder in the beaker with the solvent and cover it with a watch glass to slow down the evaporation of the solvent. Leave the beaker undisturbed and allow the solvent to flow up the paper until it reaches the 8.0 cm level.



6. The chromatogram takes about 20 to 25 minutes
 7. When the solvent in the beaker reaches the 8.0 cm level, remove the chromatogram and let it dry on a piece of paper towel without removing the staples or handling it too much with your hands.
NOTE: In case your chromatogram is still wet at the end of the laboratory session, place the chromatogram carefully into your locker for the next laboratory session.
 8. The removal of the staples and the outlining of the colored spots should only be only when the chromatogram is completely dry.
 9. Attach the chromatogram to your report form, use staples or tape.
-

Paper Chromatography Report Form

Name: _____
Date: _____
Partner: _____

1. Complete the table below:

	Number of Components	Component most attracted to station phase (color)	Component most attracted to moving phase (color)
Red Dye			
Blue Dye			
Green Dye			

2. Questions

a. What could this process be used for?

b. Name the stationary phase

c. Name the moving phase

3. Attach your chromatogram.